


Product Datasheet

Human Galectin-8 ELISA Kit (orb1728073)

Description	Human Galectin-8 ELISA Kit (96 Tests). Quantitate Human Galectin-8 in cell culture supernatants, ser...
Reactivity	Human
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-W317
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	000214
Sensitivity	25 pg/ml
Expiration Date	Please enquire.



Human
Galectin-8
PicoKine
ELISA Kit
stan...