


Product Datasheet

Human Intergrin Alpha 5/CD49E ELISA Kit (orb1728080)

Description	Human Intergrin Alpha 5/CD49E ELISA Kit (96 Tests). Quantitate Human Intergrin Alpha 5/CD49E in cell...
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F42-Y995
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	P08648
Sensitivity	75 pg/ml
Expiration Date	Please enquire.



Human
CD49E
PicoKine
ELISA Kit
standard ...