

Product Datasheet

Anti-Caspase 8/CASP8 Antibody (orb1728223)

Description

0

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IHC, IHC-Fr, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of human Caspase 8 (410-449aa VSYRNPAEGTWYIQLCQSLRERCPRGDDILTILTEVNVE), different from the related mouse and rat sequences by seven amino acids.

Form/Appearance

Lyophilized

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

55391 MW

Uniprot ID


[Q14790](#)

Dilution Range


Western blot, 0.1-0.5 µg/ml, Human, Mouse, Rat
Immunohistochemistry(Paraffin-embedded Section), 0.5-1 µg/ml, Human, Mouse, Rat, By Heat
Immunohistochemistry(Frozen Section), 0.5-1 µg/ml, Human
Immunocytochemistry, 0.5-1 µg/ml, Human
Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Flow
Cytometry
analysis of PC-
3 cells us...



Flow
Cytometry
analysis of
HeLa cells
us...



WB analysis of
Caspase8
using anti-
Caspa...