

Product Datasheet

KAP1 Antibody (orb1729745)

Description	Nuclear corepressor for KRAB domain-containing zinc finger proteins (KRAB-ZFPs). Mediates gene silen...
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide sequence around aa.827~831 (L-S-G-G-P) derived from Human TIF1 β .
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 89 kDa; Observed MW: 110 kDa
Uniprot ID	Q13263
Entrez	10155
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.