

Product Datasheet

TNF Receptor 2 Antibody (orb1729821)

Description	Receptor with high affinity for TNFSF2/TNF-alpha and approximately 5-fold lower affinity for homotri...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Immunogen	A synthetic peptide of human TNF Receptor II
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 48 kDa; Observed MW: 65 kDa
Uniprot ID	P20333
Entrez	7133
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.