

Product Datasheet

SUMO1 Antibody (orb1729860)

Description	SUMO1 Ubiquitin-like protein which can be covalently attached to target lysines as a monomer. Does n...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Sumo 1
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Uniprot ID	P63165
Entrez	7341
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.