

Product Datasheet

CCR5 Antibody (orb1729872)

Description	Receptor for a number of inflammatory CC-chemokines including CCL3/MIP-1-alpha, CCL4/MIP-1-beta and ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from human CCR5 Polyclonal
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Uniprot ID	P51681
Entrez	1234
Dilution Range	WB: 1/500-1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.