

Product Datasheet

MAP4K2/3 Antibody (orb1729876)

Description	Serine/threonine-protein kinase which acts as an essential component of the MAP kinase signal transd...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from human GCK/GLK. at AA range: 221-270
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 92 kDa; Observed MW: 100 kDa
Uniprot ID	Q8IVH8 , Q12851
Entrez	8491 , 5871
Dilution Range	WB: 1/500-1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.