

Product Datasheet

PSPN Antibody (orb1729928)

Description	The protein encoded by this gene is a neurotrophic factor, belonging to the GDNF family. Neurotrophi...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human PSPN. AA range:101-150
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Uniprot ID	O60542
Entrez	5623
Dilution Range	WB: 1/500-1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.