

Product Datasheet

Histone H1 Antibody (orb1729999)

Description	Histone H1 protein binds to linker DNA between nucleosomes forming the macromolecular structure know...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Histone H1. AA range:1-50
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 22 kDa; Observed MW: 30 kDa
Uniprot ID	P16401, P16402, P10412
Entrez	3007, 3008, 3009
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; ELISA: 1/10000
Expiration Date	12 months from date of receipt.