

Product Datasheet

Mannose Receptor Antibody (orb1730013)

Description	Mediates the endocytosis of glycoproteins by macrophages. Binds both sulfated and non-sulfated polys...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human MRC1. AA range:341-390
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 166 kDa; Observed MW: 170 kDa
Uniprot ID	P22897
Entrez	4360
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.