

## Product Datasheet

### CRTH2 Antibody (orb1730021)

|                            |                                                                                                                           |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Receptor for prostaglandin D2 (PGD2).                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                    |
| <b>Reactivity</b>          | Human, Mouse                                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                                              |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                                            |
| <b>Immunogen</b>           | The antiserum was produced against synthesized peptide derived from the Internal region of human PTGDR2. AA range:161-210 |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.                                      |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                  |
| <b>Note</b>                | For research use only                                                                                                     |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                                |
| <b>MW</b>                  | Calculated MW: 43 kDa; Observed MW: 43 kDa                                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q9Y5Y4</a>                                                                                                    |
| <b>Entrez</b>              | <a href="#">11251</a>                                                                                                     |
| <b>Dilution Range</b>      | WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                           |