

Product Datasheet

CXCR2 Antibody (orb1730022)

Description	The protein encoded by this gene is a member of the G-protein-coupled receptor family. This protein ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CXCR2. AA range:1-50
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Uniprot ID	P25025
Entrez	3579
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.