

Product Datasheet

EYA1/4 Antibody (orb1730033)

Description	Functions both as protein phosphatase and as transcriptional coactivator for SIX1, and probably also...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide from human protein at AA range: 271-320
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	095677, Q99502
Entrez	2070, 2138
Dilution Range	IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.