

Product Datasheet

MYH Antibody (orb1730066)

Description	Muscle contraction.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MYH-pan around the non-acetylation site of Lys1394. AA range:1351-1400
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 223 kDa; Observed MW: 223 kDa
Uniprot ID	Q9UKX2 , Q9Y623 , P12882 , P11055 , P13533 , P13535 , P12883
Entrez	4619 , 4621 , 4626 , 4624 , 4620 , 4625 , 4622
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.