

Product Datasheet

IGF1 Receptor Antibody (orb1730079)

Description	This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Peptide sequence around aa.1159~1163 (D-I-Y-E-T) derived from Human IGF-1R .
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 155 kDa; Observed MW: 95,200 kDa
Uniprot ID	P08069
Entrez	3480
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.