

Product Datasheet

NT5DC3 Antibody (orb1730112)

Description	Predicted to enable 5'-nucleotidase activity. Predicted to be involved in dephosphorylation. Part of...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized peptide derived from human NT5DC3
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 63 kDa; Observed MW: 60 kDa
Uniprot ID	Q86UY8
Entrez	51559
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.