

Product Datasheet

SGLT1 Antibody (orb1730121)

Description	Actively transports glucose into cells by Na ⁺ cotransport with a Na ⁺ to glucose coupling ratio of 2:...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SGLT-1. AA range:525-574
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Uniprot ID	P13866
Entrez	6523
Dilution Range	WB: 1/500-1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.