

Product Datasheet

Myomesin 2 Antibody (orb1730136)

Description	Major component of the vertebrate myofibrillar M band. Binds myosin, titin, and light meromyosin. Th...
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MYOM2. AA range:612-661
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 165 kDa; Observed MW: 165 kDa
Uniprot ID	P54296
Entrez	9172
Dilution Range	WB: 1/500-1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.