

Product Datasheet

Metalloproteinase Inhibitor 1 Antibody (orb1730156)

Description

TIMP1 Complexes with

collagenases) and irreversibly inactivates them by ...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IHC, WB

Applications
Immunogen

The antiserum was produced against synthesized peptide derived from human TIMP1. AA range:61-110

Preservatives

Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

Calculated MW: 23 kDa; Observed MW: 23 kDa

Uniprot ID

P01033

Entrez

7076

Dilution Range

WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; ELISA: 1/10000

Expiration Date

12 months from date of receipt.