

Product Datasheet

KA1 Antibody (orb1730191)

Description	Receptor for glutamate. L-glutamate acts as an excitatory neurotransmitter at many synapses in the c...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human KA1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 107 kDa; Observed MW: 107 kDa
Uniprot ID	Q16099
Entrez	2900
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.