

Product Datasheet

SLC9A1 Antibody (orb1730217)

Description	Involved in pH regulation to eliminate acids generated by active metabolism or to counter adverse en...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human Sodium / Hydrogen Exchanger 1
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 91 kDa; Observed MW: 130 kDa
Uniprot ID	P19634
Entrez	6548
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.