

Product Datasheet

IGFBP7 Antibody (orb1730236)

Description	Binds IGF-I and IGF-II with a relatively low affinity. Stimulates prostacyclin (PGI ₂) production. St...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, WB
Immunogen	A synthesized peptide derived from human IGFBP7
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Uniprot ID	Q16270
Entrez	3490
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.