

Product Datasheet

STING Antibody (orb1730266)

Description	Facilitator of innate immune signaling that promotes the production of type I interferon (IFN-alpha ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	A synthesized peptide derived from human TMEM173
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 42 kDa; Observed MW: 37 kDa
Uniprot ID	Q86WV6
Entrez	340061
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.