

Product Datasheet

Substance P Antibody (orb1730387)

Description	Tachykinins are active peptides which excite neurons, evoke behavioral responses, are potent vasodil...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthesized peptide derived from human Substance P
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	P20366
Entrez	6863
Dilution Range	IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.