

Product Datasheet

Cyclophilin B Antibody (orb1730517)

Description	PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imid...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human Cyclophilin B
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 24 kDa; Observed MW: 18 kDa
Uniprot ID	P23284
Entrez	5479
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.