

Product Datasheet

SCARB1 Antibody (orb1730522)

Description	Receptor for different ligands such as phospholipids, cholesterol ester, lipoproteins, phosphatidyls...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, WB
Immunogen	A synthesized peptide derived from human Scavenging Receptor SR-BI
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 61 kDa; Observed MW: 80 kDa
Uniprot ID	Q8WTV0
Entrez	949
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.