

Product Datasheet

Rad51 Antibody (orb1730543)

Description Rad51 participates in a common DNA damage response pathway associated with the activation of homolog...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IHC, IP, WB

Immunogen A synthesized peptide derived from human Rad51

Preservatives Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW Calculated MW: 37 kDa; Observed MW: 37 kDa

Uniprot ID [Q06609](#)

Entrez [5888](#)

Dilution Range IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100

Expiration Date 12 months from date of receipt.