

Product Datasheet

COL4A1 Antibody (orb1730555)

Description	This gene encodes a type IV coll...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from COL4A1. at AA range: 1428-1443
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	161kDa
Entrez	1282
Dilution Range	WB: 1/500-1/2000; ELISA: 1/10000-1/20000
Expiration Date	12 months from date of receipt.