

Product Datasheet

Acetyl-Lysine Antibody (orb1730558)

Description	Acetylation of lysine residues within proteins has emerged as an important mechanism used by cells t...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein of Acetyl Kine
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Refer to figures
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.