

Product Datasheet

Phospho-JNK (Thr183/Tyr185) Antibody (orb1730565)

Description

The protein encoded by this gene is a

kinases act as an integra...

Species/Host

Rabbit

Reactivity

Gallus, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IHC, WB

Applications
Immunogen

The antiserum was produced against synthesized peptide derived from human JNK1/2/3 around the phosphorylation site of Thr183 and Tyr185. AA range:151-200

Preservatives

Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

Calculated MW: 48 kDa; Observed MW: 46,54 kDa

Uniprot ID

P45983

Entrez

5599

Dilution Range

WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; ELISA: 1/10000

Expiration Date

12 months from date of receipt.