

Product Datasheet

Phospho-Smad2/3 (Thr8) Antibody (orb1730576)

Description	The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene produ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Smad2/3 around the phosphorylation site of Thr8. AA range:1-50
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 52 kDa; Observed MW: 48-52~60 kDa
Uniprot ID	P84022, Q15796
Entrez	4087, 4088
Dilution Range	WB: 1/500-1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.