

Product Datasheet

Phospho-STAT3 (Ser727) Antibody (orb1730578)

Description	The protein encoded by this gene is a member of the STAT protein family. In response to cytokines an...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Phospho-STAT3 (S727)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 88 kDa; Observed MW: 98 kDa
Uniprot ID	P40763
Entrez	6774
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.