

Product Datasheet

Phospho-p38 (Tyr323) Antibody (orb1730595)

Description	Responds to activation by environmental stress, pro-inflammatory cytokines and lipopolysaccharide (L...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human MAPK14
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 41 kDa; Observed MW: 35 kDa
Uniprot ID	Q16539
Entrez	1432
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.