

Product Datasheet

Phospho-FRS2 (Tyr436) Antibody (orb1730611)

Description	Adapter protein that links FGR and NGF receptors to downstream signaling pathways. Involved in the a...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide of human FRS2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 57 kDa; Observed MW: 65 kDa
Uniprot ID	Q8WU20
Entrez	10818
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.