

Product Datasheet

Phospho-STAT6 (Tyr641) Antibody (orb1730623)

Description	STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological respo...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide of human STAT6
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 94 kDa; Observed MW: 110 kDa
Uniprot ID	P42226
Entrez	6778
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.