

Product Datasheet

Phospho-NMDAR2A/2B (Tyr1246/Tyr1252) Antibody (orb1730629)

Description	NMDA receptor subtype of glutamate-gated ion channels possesses high calcium permeability and voltag...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Immunogen	The antiserum was produced against synthesized peptide derived from human NMDAR2A/B around the phosphorylation site of Tyr1246/1252. AA range:1216-1265
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q13224, Q12879
Entrez	2904, 2903
Dilution Range	IHC: 1/50-1/100; ICC: 1/50-1/200; ELISA: 1/10000
Expiration Date	12 months from date of receipt.