

Product Datasheet

**Phospho-KCNQ2/3/4/5
(Thr217/Thr246/Thr223/Thr251) Antibody
(orb1730631)**

Description

The M channel is a slowly

potassium channel that plays a critical role i...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Immunogen

The antiserum was produced against synthesized peptide derived from human Kv7.3/KCNQ3 around the phosphorylation site of Thr246. AA range:191-240

Preservatives

Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

Uniprot ID

P56696, O43526, O43525, Q9NR82

Entrez

3785, 9132, 3786, 56479

Dilution Range

IHC: 1/50-1/100; ELISA: 1/10000

Expiration Date

12 months from date of receipt.