

Product Datasheet

Phospho-p107 (Thr369) Antibody (orb1730635)

Description	Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintai...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	The antiserum was produced against synthesized peptide derived from human RBL1 around the phosphorylation site of Thr369. AA range:335-384
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	P28749
Entrez	5933
Dilution Range	IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.