

Product Datasheet

Phospho-JAK3 (Tyr785) Antibody (orb1730647)

Description	Tyrosine kinase of the non-receptor type, involved in the interleukin-2 and interleukin-4 signaling ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human JAK3 around the phosphorylation site of Tyr785. AA range:751-800
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 125 kDa; Observed MW: 125 kDa
Uniprot ID	P52333
Entrez	3718
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.