

Product Datasheet

Phospho-STAT3 (Ser727) Antibody (orb1730649)

Description	The protein encoded by this gene is a member of the STAT protein family. In response to cytokines an...
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, IP, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human STAT3 around the phosphorylation site of Ser727. AA range:694-743
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 88 kDa; Observed MW: 85 kDa
Uniprot ID	P40763
Entrez	6774
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; ELISA: 1/10000
Expiration Date	12 months from date of receipt.