

Product Datasheet

Phospho-Src (Tyr419) Antibody (orb1730689)

Description	Src activity is regulated by tyrosine phosphorylation at two sites, but with opposing effects. While...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthesized peptide derived from human Phospho-Src (Y419)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Uniprot ID	P12931
Entrez	6714
Dilution Range	WB: 1/500-1/1000; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.