

Product Datasheet

Phospho-eIF4E (Ser209) Antibody (orb1730691)

Description	eIF4E, a protein modulates translation of maternal mRNAs in early embryos before the onset of zygoti...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Phospho-eIF4E (S209)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Uniprot ID	P06730
Entrez	1977
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.