

Product Datasheet

Phospho-RSK1 p90 (Thr359/Ser363) Antibody (orb1730694)

Description	Rsk1 is a member of a family of 90kDa ribosomal protein S6 kinases, which includes Rsk1, Rsk2 and Rs...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Phospho-RSK1 p90 (T359+S363)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 83 kDa; Observed MW: 90 kDa
Uniprot ID	Q15418
Entrez	6195
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.