

Product Datasheet

Acetyl-Histone H4 (Lys16) Antibody (orb1730723)

Description	The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the D...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic Peptide of Histone H4 (Acetyl Lys16)
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 11 kDa; Observed MW: 14 kDa
Uniprot ID	P62805
Entrez	121504
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.