

Product Datasheet

DAP kinase 3 Antibody (orb1730731)

Description	Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpres...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human ZIP Kinase
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Uniprot ID	043293
Entrez	1613
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.