

Product Datasheet

Transcription Initiation Factor IIB Antibody (orb1730765)

Description	General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA poly...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human TFIIIB
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Uniprot ID	Q00403
Entrez	2959
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.