

Product Datasheet

NUP98 Antibody (orb1730870)

Description	Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Nup98 and Nup96 are invo...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human NUP98
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 198 kDa; Observed MW: 98 kDa
Uniprot ID	P52948
Entrez	4928
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.