

Product Datasheet

XPD Antibody (orb1731185)

Description	ATP-dependent 5'-3' DNA helicase, component of the core-TFIIF basal transcription factor. Involved i...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	A synthetic peptide of human XPD
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 87 kDa; Observed MW: 80 kDa
Uniprot ID	P18074
Entrez	2068
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.