

Product Datasheet

TPPP Antibody (orb1731232)

Description	May play a role in the polymerization of tubulin into microtubules, microtubule bundling and the sta...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	A synthetic peptide of human Tubulin Polymerization Promoting Protein
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Uniprot ID	O94811
Entrez	11076
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.